



AUSTROMAT™ 220



The DEKEMA for starters. No compromise on results.

compact - precise - full-featured

- M Stand-alone furnace with optional vacuum pump
- 13 kg stable
- In long lasting quality, using noble, light materials
- ➢ Firing, glazing, crystallization and infiltration with and without vacuum
- Place power perform

DEKEMA. Chamber.

- Well known from the high quality AUSTROMAT™ Series 6 Firing Furnace
- Designed decades ago, operated for decades, optimized over decades
- Unique, homogeneous energy distribution
- ▶ Fiber insulation of purest quality in highest precision design







DEKEMA. Heating.

- Well known from the high quality
 AUSTROMAT™ Series 6 Firing Furnace
- Manufactured in house by DEKEMA
- Unsurpassed durability and longevity

DEKEMA. Firing table.

- Well known from the high quality
 AUSTROMAT™ Series 6 Firing Furnace
- Physically defined fiber insulation table of outstanding purity
- Usable with common firing trays

DEKEMA. Casing.

- ➢ Fan free, integrated convection cooling
- ⚠ Corrosion free, lightweight aluminum housing
- ✗ Kept in matt black



DEKEMA. Thermocouple & Control.

- Well known from the high quality
 AUSTROMAT™ Series 6 Firing Furnace
- ▶ Platinum/platinum-rhodium thermocouple of the purest quality and a controlled accuracy far below 1°C
- ⚠ Temperature calibration at only one measuring point through the DEKEMA calibration set, unmatched in precision

DEKEMA. Lift system.

- ▲ Smooth lift movement
- ▲ Durable stainless steel design
- Freely selectable drying position

DEKEMA. User Interface.

- ▲ Language independent icon surface
- ✗ Touch operation with 100 programs
- Independent vacuum start







Optional. For maximum flexibility.

DEKEMA. Vacuum pump.

- Low contamination Trenn Coat™ coating
- > Straight air flow
- Low maintenance design
- ✗ Specific shell



DEKEMA. Device equipment.

	AUSTROMAT™ 220	AUSTROMAT™ 624	AUSTROMAT™ 624 <i>i</i> ™
Size (WxHxD)	10" x 21" x 10"	15"x26"x14"	15"x26"x14"
Weight	28,5 lbs	38,5 lbs	42 lbs
Power Supply	230 V/50-60 Hz	230 V/50-60 Hz	230 V/50 – 60 Hz
Maximum Power Supply	1.500 W	1.500 W	1.500 W
T_{max}	1.100 °C	1.200 °C	1.200 °C
Touchscreen	4,3 Zoll, color	7,5 Zoll, color	10,4 Zoll, color
USB interface / Ethernet interface	-/-	/ optional	<i>)</i>),
DEKEMA Autodry™	-	<i>I</i> II	<i>)</i>),
Housing: Aluminium	<i>I</i> II	<i>I</i> II	<i>)</i>),
Corporate Design Customization	-	optional	optional
Quartz glass heating element: DEKEMA specific heating element coil guided in pure quartz glass	<i>)</i>	I II	<i>)</i>),
Fiber insulation: DEKEMA specific fiber insulation	<i>I</i> II	<i>I</i> II	<i>)</i>),
Fiber insulation table: For firing porcelain, adapted to the dimentions of the furnace)h	<i>)</i>	<i>)</i>),
DEKEMA Vakuum pump	optional	<i>)</i> }	<i>)</i>),





DEKEMA. Software Features.

	AUSTROMAT™ 220	AUSTROMAT™ 624	AUSTROMAT™ 624 <i>i</i> ™
DEKEMA Internet database access for registered users		I h	<i>)</i>),
DEKEMA Internet database furnace connectivity	-	optional	<i>)</i>),
Internet program streaming	-	-	<i>)</i>),
Cloud: Data management in the DEKEMA Circle	-	-	<i>)</i>),
IoT (Internet of Things: If the DEKEMA furnace is registered in the DEKEMA online portal, the user benefits from a secure backup of the process data and a fast service to minimize downtime.)	-	-	<i>)</i>),
iDREAM: VNC-Viewer for remote control	-	optional	<i>)</i>),
X-DREAM: Computer program for querying the furnace state and emailing QM data	-	optional	<i>)</i>),
QM (quality management): Record of process files	-	optional	
USB: Transferring firing programs (dds)	-))	<i>)</i>),
USB: Transferring fotos (jpeg)	-		
USB: Storage of QM files (dfl)	-	optional	<i>)</i>),
USB: Storage of Printouts	-		
USB: Backup/Update	-	I II	<i>)</i>),
Internet: Backup / Update	-	optional	
FTP-Server	-	optional	<i>)</i>),
FTP-Client	-	optional	<i>)</i>),
OPC: Software server to incorporate quality management into ERP systems	-	optional	optional
Programming languages	-	Classic, Easy, Professional	Classic, Easy, Professional, Professional Plus
Language-independent program management	<i>)</i>),	-	-







Der Ofen.™

DEKEMA Dental-Keramiköfen GmbH

Industriestraße 22 • 83395 Freilassing • Germany ©+49-8654-46390 **(a)** +49-8654-66195 info@dekema.com • www.dekema.com





